

Product Datasheet

Anti-Epstein-Barr Virus BcLF1 Antibody (orb2282882)

Description	Anti-BcLF1 Antibody catalog # orb2282882. This antibody reacts with Epstein-Barr Virus.
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.