

Product Datasheet

Mouse Anti-Rabbit IgG mAb/Gold (orb2283905)

Description	Mouse Anti-Rabbit IgG mAb/Gold
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Gold
Immunogen	Native Rabbit IgG full length
Preservatives	20mM TBS (pH8.0) with 1% rAlbumin and 0.02% Proclin300.
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.