

Product Datasheet

phospho-TSC1 (Ser505)/BF555 Antibody (orb2283951)

Description	phospho-TSC1 (Ser505)/BF555 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	BF555
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TSC1 around the phosphorylation site of Ser505 FD(p-S)PF
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92574
Expiration Date	12 months from date of receipt.