

Product Datasheet

IL2RA/CD25/PerCP-Cy5.5 Antibody (orb2284240)

Description	IL2RA/CD25/PerCP-Cy5.5 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human IL2RA/CD25
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P01590
Expiration Date	12 months from date of receipt.