

Product Datasheet

SARS-CoV-2 Antibody (Biotin) (orb2284273)

Description	SARS-CoV-2 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Spike S1 Protein (14-685/1213aa)
Target	SARS-CoV-2
Preservatives	0.01M PBS(pH7.4)
Concentration	Lot Dependent
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.