

Product Datasheet

CDC25A Antibody (BF488) (orb2284293)

Description	CDC25A Antibody (BF488)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	BF488
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human cdc25A around the phosphorylation site of Ser76 MG(p-S)SE
Target	CDC25A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P30304
Expiration Date	12 months from date of receipt.