

Product Datasheet

RPS6 Antibody (BF350) (orb2284445)

| | |
|----------------------------|---|
| Description | RPS6 Antibody (BF350) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | BF350 |
| Tested Applications | FC |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human RPS6 around the phosphorylation site of Ser240/244 RA(p-S)TSK(p-S)ES |
| Target | RPS6 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P62753 |
| Expiration Date | 12 months from date of receipt. |