

Product Datasheet

RPS6 Antibody (BF488) (orb2284446)

Description	RPS6 Antibody (BF488)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF488
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RPS6 around the phosphorylation site of Ser240/244 RA(p-S)TSK(p-S)ES
Target	RPS6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62753
Expiration Date	12 months from date of receipt.