

## Product Datasheet

# MAP1LC3C (Human) Recombinant Protein (P01) (orb2284557)

<b>Description</b>	Human MAP1LC3C full-length ORF ( NP_001004343.1, 1 a.a. - 147 a.a.) recombinant protein with GST-tag...
<b>Tested Applications</b>	AP, Array, ELISA, WB
<b>Preservatives</b>	50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer.
<b>Storage</b>	Store at -80°C. Aliquot to avoid repeated freezing and thawing.
<b>Tag</b>	GST
<b>Note</b>	For research use only
<b>Application notes</b>	Best use within three months from the date of receipt of this protein.
<b>Protein Sequence</b>	MPPPQKIPSVRPFKQRKSLAIRQEEVAGIRAKFPNKIPVVVERYPRETFLPPLDKTK FLVPQELTMTQFLSIIRSRMVLRATEAFYLLVNNKSLVMSATMAEIYRDYKDEDG FVYMTYASQETFGCLESAAAPRDGSSLEDRPCNPL
<b>Clonality</b>	Recombinant
<b>NCBI</b>	<a href="https://www.ncbi.nlm.nih.gov/nuccore/NP_001004343.1">NP_001004343.1</a>
<b>Expiration Date</b>	6 months from date of receipt.