

## Product Datasheet

# YTHDF3 (Human) Recombinant Protein (P01) (orb2284617)

<b>Description</b>	Human YTHDF3 full-length ORF ( NP_689971.4, 1 a.a. - 585 a.a.) recombinant protein with GST-tag at N...
<b>Tested Applications</b>	AP, Array, ELISA, WB
<b>Preservatives</b>	50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer.
<b>Storage</b>	Store at -80°C. Aliquot to avoid repeated freezing and thawing.
<b>Tag</b>	GST
<b>Note</b>	For research use only
<b>Application notes</b>	Best use within three months from the date of receipt of this protein.
<b>Protein Sequence</b>	MSATSVDQRPKGQGKNSVQNGSIHQKDAVNDDDFEPYLSSQTNQSNSYPPM SDPYMPSYYAPSIGFPYSLGEAAWSTAGDQMPYLTTYGQMSNGEHYIPDGVF SQPGALGNTPPFLGQHGFNFFPGNADFSTWGTSGSQGQSTQSSAYSSSYGYPP SSLGRAITDQGAGFGNDTSLKVPGISSIEQGMTGLKIGGDLTAAVTKTVGTALSS SGMTSIATNSVPPVSSAAPKPTSWAAIARKPAKQPQKPKGNVIGGSAVPPPI KHNMNIGTWDEKGSVVKAPPTQPVLPQTIIQQPQLIQQPPPLVQSQLPQQQPQ PQPQQQGPQPQAQPHQVQPQQQLQNRWVAPRNRGAGFNQNGAGSENF GLGVVPSASPSSEVHPVLEKLNKAINNYPKDFDWNLNKGRVFIKSYSEDDIHR SIKYSIWCSTEHGKRLDAAAYRSLNGKGPLYLLFSVNGSGHFCGVAEMKSVVDY NAYAGVWSQDKWKGKFEVKWIFVKDVPNNQLRHIRENNDNKPVNTSRDTQE VPLEKAKQVLKIIATFKHTTSIFDDFAHYEKRQEEEEAMRRERNRNKQ
<b>Clonality</b>	Recombinant
<b>NCBI</b>	<a href="#">NP_689971.4</a>
<b>Expiration Date</b>	6 months from date of receipt.