

Product Datasheet

RIOK1 (Human) Recombinant Protein (P01) (orb2284899)

Description	Human RIOK1 full-length ORF (NP_113668.2, 1 a.a. - 568 a.a.) recombinant protein with GST-tag at N-...
Tested Applications	AP, Array, ELISA, WB
Preservatives	50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer.
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.
Tag	GST
Note	For research use only
Application notes	Best use within three months from the date of receipt of this protein.
Protein Sequence	MDYRLLMSRVVPGQFDDADSSDSENRLKTVKEKDDILFEDLQDNVNENGE EIEDEEEEGYDDDDDDWDWDEGVGKLAKGYVWNGGSPQANRQTSDDSSAK MSTPADKVLRFENKINLDKLNVTDSVINKVTEKSRQKEADMYRIKDKADRATVE QVLDPRTRMILFKMLTRGIITEINGCISTGKEANVYHASTANGESRAIKIYKTSILVF KDRDKYVSGEFRFRHGYCKGNPRKMKVKTWAEKEMRNLRNLNTAEIPCPEPIMLR SHVLVMSFIGKDDMPAPLLKNVQLSESKARELYLQVIQYMRRMYQDARLVHADL SEFNMLYHGGGVYIIDVSQSVEHDHPHALEFLRKDCANVNDFFMRHRSVAVMTV RELFVFTDPSITHENMDAYLSKAMEIASQRTKEERSSQDHVDEEVFKRAYIPRTL NEVKNYERDMDIIMKLKEEDMAMNAQQDNILYQTVTGLKDLGSGVQKVPALLEN QVEERTCSDSEDIGSSECSDTDSEEQGDHARPKHHTTDPDIDKKERKKMVKEAQ REKRKNKIPKHVKKRKEKTAKTKKGK
Clonality	Recombinant
NCBI	NP_113668.2
Expiration Date	6 months from date of receipt.