

Product Datasheet

VPS33B (Human) Recombinant Protein (P01) (orb2285480)

Description	Human VPS33B full-length ORF (AAH16445.1, 1 a.a. - 617 a.a.) recombinant protein with GST-tag at N-...
Tested Applications	AP, Array, ELISA, WB
Preservatives	50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer.
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.
Tag	GST
Note	For research use only
Application notes	Best use within three months from the date of receipt of this protein.
Protein Sequence	MAFPHRPDAPELPDFSMKRLARDQLIYLLEQLPGKKDLFIEADLMSPLDRIVNS ILKQHEVDKLYKVENKPALSSNEQLCFLVRPRIKNMRYIASLVNADKLAGRTRKYK VIFSPQKFYACEMVLEEEGIYGDVSCDEWAFSLLPLDLDLLSMELPEFFRDYFLEG DQRWINTVAQALHLLSTLYGPFPCYIGRCAKMAYELWRNLEEEEDGETKGRR PEIGHIFLLDRDVFVTALCSQVVYEGLVDDTFRIKCGSVDFGPEVTSSDKSLKVL LNAEDKVFNEIRNEHFSNVFGFLSQKARNLQAQYDRRRGMDIKQMKNFVSQELK GLKQEHRLLSLHIGACESIMKKKTKQDFQELIKTEHALLEGFNIRESTSYIEEHIDR QVSPIESLRLMCLLSITENGLIPKDYRSLKTQYLQSYGPEHLLTFSNLRRAGLLTEQ APGDTLTAVESKVSCLVTDKAAGKITDAFSSLAKRSNFRAISKLNLIIPRVDGEYD LKVPRDMAYVFGAYVPLSRIIEQVLERRSWQGLDEVVRLNCSDFAFDTMTK EDKASSESLRLLVFLGGCTFSEISALRFLGREKGYRFIFLTTAVTNSARLMEAMS EVKA
Clonality	Recombinant
NCBI	AAH16445.1
Expiration Date	6 months from date of receipt.