

Product Datasheet

VEGFB (Human) Recombinant Protein (P01) (orb2286323)

Description	Human VEGFB full-length ORF (AAH08818.1, 1 a.a. - 207 a.a.) recombinant protein with GST-tag at N-t...
Tested Applications	AP, Array, ELISA, WB
Preservatives	50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer.
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.
Tag	GST
Note	For research use only
Application notes	Best use within three months from the date of receipt of this protein.
Protein Sequence	MSPLLRLLLLLAALLQLAPAQAPVSQPDAPGHQRKVVSWIDVYTRATCQPREVVV PLTVELMGTVAKQLVPSCVTVQRCGGCCPDDGLECVPTGQHQVRMQILMIRYP SSQLGEMSLEEHSQCECRPKKKDSAVKPDRAATPHHRPQPRSVPGWDSAPGAP SPADITHPTPAPGPSAHAAPSTTSALTPGPAAAAADAAASSVAKGGA
Clonality	Recombinant
NCBI	AAH08818.1
Expiration Date	6 months from date of receipt.