

## Product Datasheet

# MAP1LC3C MaxPab mouse polyclonal antibody (B01) (orb2290555)

<b>Description</b>	Mouse polyclonal antibody raised against a full-length human MAP1LC3C protein.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	MAP1LC3C (NP_001004343, 1 a.a. ~ 147 a.a) full-length human protein.
<b>Preservatives</b>	No additive
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	For IHC and IF applications, antibody purification with Protein A will be needed prior to use.
<b>Protein Sequence</b>	MPPPQKIPSVRPFKQRKSLAIRQEEVAGIRAKFPNKIPVVVERYPRETFLPPLDKTK FLVPQELTMTQFLSIIRSRMVLRLATEAFYLLVNNKSLVMSATMAEIYRDYKDEDG FVYMTYASQETFGCLESAAPRDGSSLEDRPCNPL
<b>Clonality</b>	Polyclonal
<b>NCBI</b>	<a href="#">NP_001004343</a>
<b>Expiration Date</b>	12 months from date of receipt.