

---

## Product Datasheet

**ADPGK monoclonal antibody (M01), clone 1E4  
(orb2290776)**

**Description**

Mouse monoclonal antibody raised against a partial recombinant


**Species/Host**

Mouse

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, WB

**Immunogen**

ADPGK (NP\_112574, 425 a.a. ~ 496 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Preservatives**

In 1x PBS, pH 7.4

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Protein Sequence**

 SLRAPQEFMTSHSEAGSRIVLNPKNKPVVEWHREGISFHFTPVLVCKDPIRTVGLG  
DAISAEGLFYSEVHPHY

**Isotype**


IgG2a Kappa

**Clonality**


Monoclonal

**NCBI**
[NP\\_112574](https://www.ncbi.nlm.nih.gov/nuccore/NP_112574)
**Expiration Date**


12 months from date of receipt.



ADPGK monoclonal antibody (M01), clone 1...



ADPGK monoclonal antibody (M01), clone 1...



ADPGK monoclonal antibody (M01), clone 1...