


Product Datasheet

Histone H3K79me3 antibody (orb229082)

Description	Rabbit polyclonal antibody to Histone H3K79me3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC-P, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:20-1:50, IP: 1:20-1:50, ChIP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell line ...