

Product Datasheet

NXT1 purified MaxPab rabbit polyclonal antibody (D01P) (orb2291213)

Description	Rabbit polyclonal antibody raised against a full-length human NXT1 protein.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	NXT1 (AAH00759.1, 1 a.a. ~ 140 a.a) full-length human protein.
Preservatives	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MASVDFKTYVDQACRAAEFVNYYTTMDKRRRLLSRLYMGATLVWNGNAVS GQESSSEFFEMLPSEFQISVVDCQPVHDEATPSQTTVLVVICGSVKFEGNKQRD FNQNFILTAQASPSNTVWKIASDCFRFQDWAS
Clonality	Polyclonal
NCBI	AAH00759.1
Expiration Date	12 months from date of receipt.