



Product Datasheet

HEY1 monoclonal antibody (M02), clone 2F10 (orb2291330)


Description	Mouse monoclonal antibody raised against a partial recombinant HEY1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	HEY1 (NP_036390, 121 a.a. ~ 220 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Preservatives	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	DYRSLGFRECLAEVARYLSIIIEGLDASDPLRVRLVSHLNYYASQREAAASGAHAGL GHIPWGTVFGHHPHIAHPLLLPQNGHGNAGTTASPTEPHHQGRGLG
Isotype	IgG2a Kappa
Clonality	Monoclonal
NCBI	NP_036390
Expiration Date	12 months from date of receipt.



Detection limit for recombinant...
GST tagg...



Western Blot analysis of HEY1 expression...



Western Blot detection against Immunogen...