



Product Datasheet

SAE1 monoclonal antibody (M01), clone 1G4-1G5 (orb2291630)


Description	Mouse monoclonal antibody raised against a full length recombinant SAE1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	SAE1 (AAH00344, 1 a.a. ~ 346 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Preservatives	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MVEKEEAGGGISEEEAAQYDRQIRLWGLEAQKRLRASRVLLVGLKGLGAEIAKNL ILAGVKGLTMLDHEQVTPEDPGAQFLIRTGSVGRNRAEASLERAQNLNPMVDVK VDTEIEKKPESFFTQFDAVCLTCCSRDVIVKVDQICHKNSIKFFTGDVFGYHGYT FANLGEHEFVEEKTQVAKVSQGVEDGPDTKRAKLDSSSETTMVKKKVVFCPVKEA LEVDWSSEKAKAALKRTTSDYFLLQVLLKFRDVKGRDPSSDYEEDSELLQIRN DVLDSLGISPDLLPEDFVRYCFSEMAPVCAVVGILAQEIVKALSQRDPPHNNFF FDGMKGNNGIVECLGPK
Isotype	IgG1 lambda
Clonality	Monoclonal
NCBI	AAH00344
Expiration Date	12 months from date of receipt.



Western Blot detection against Immunogen...



Detection limit for recombinant GST tagg...



Immunofluorescence of monoclonal antibody...