



## Product Datasheet

# SULT1A1 monoclonal antibody (M01A), clone 1F8 (orb2292110)


<b>Description</b>	Mouse monoclonal antibody raised against a full length recombinant SULT1A1.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IP, WB
<b>Immunogen</b>	SULT1A1 (AAH00923, 1 a.a. ~ 295 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>Preservatives</b>	In ascites fluid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MELIQDTSRPPLEYVKGVPLIKYFAEALGPLQSFQARPDDLSTYPKSGTTWVSQI LDMIIYQGGDLEKCHRAPIFMRVPFLEFKAPGIPSGMETLKDTPAPRLLKTHLPLALL PQTLLDQKVKVYVARNAKDVAVSYYHFYHMAKVHPEPGTWDSFLEKFMVGEV SYGSWYQHVVQEWELSRTHPVLYLFYEDMKENPKREIQKILEFVGHSLPEETVD FVVQHTSFKEMKKNPMTNYTTVPQEFMDHSISPFMRKGMAGDWKTTFTVAQNE RFDADYAEKMAGCSLSFRSEL
<b>Isotype</b>	IgM kappa
<b>Clonality</b>	Monoclonal
<b>NCBI</b>	<a href="#">AAH00923</a>
<b>Expiration Date</b>	12 months from date of receipt.



Immunoprecipitation of SULT1A1 transfect...



Western Blot analysis of SULT1A1 express...



Western Blot detection against Immunogen...