



## Product Datasheet

### SCARB2 polyclonal antibody (A01) (orb2294808)


<b>Description</b>	Mouse polyclonal antibody raised against a partial recombinant SCARB2.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	SCARB2 (NP_005497, 339 a.a. ~ 437 a.a) partial recombinant protein with GST tag.
<b>Preservatives</b>	50 % glycerol
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	FPHFYQADERFVSAIEGMHPNQEDHETFVDINPLTGILKAAKRFQINIYVKKLDDF VETGDIRTMVFPVMYL NESVHIDKETASRLKSMINTTLITN
<b>Clonality</b>	Polyclonal
<b>NCBI</b>	<a href="#">NP_005497</a>
<b>Expiration Date</b>	12 months from date of receipt.



SCARB2 polyclonal antibody (A01).  
Wester...



Western Blot analysis of SCARB2 expressi...



Western Blot detection against Immunogen...