



Product Datasheet

B2M monoclonal antibody (M01J), clone 3F9-2C2 (orb2295607)


Description	Mouse monoclonal antibody raised against a full-length recombinant B2M. This product is belong to Cel...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, PLA, WB
Immunogen	B2M (AAH32589.1, 1 a.a. ~ 119 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Preservatives	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MSRSVALAVLALLSLSGLEAIQRTPKIQVYSRHPAENGKSNFLNCYVSGFHPSDIEVDLLKNGERIEKVEHSDLSFSKDWFSYLLYYTEFTPTTEKDEYACRVNHVTLSPKIVKWDRDM
Isotype	IgG2b Kappa
Clonality	Monoclonal
NCBI	AAH32589.1
Expiration Date	12 months from date of receipt.



B2M monoclonal antibody (M01J), clone 3F...



Detection limit for recombinant GST tagg...



Immunofluorescence of monoclonal antibod...