

Product Datasheet

ARG2 purified MaxPab mouse polyclonal antibody (B01P) (orb2295888)


Description	Mouse polyclonal antibody raised against a full-length human ARG2 protein.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	ARG2 (AAH29050.1, 1 a.a. ~ 354 a.a) full-length human protein.
Preservatives	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MSLRGSLSRLLQTRVHSILKKS VHSVAVIGAPFSQGQKRKGVEHGPA AIREAGLMKRLSSLG CHLKDFGDL SFTVPVKDDL YNNLIVNPRSVGLANQELAEVVSRAVSDGYSCVTLGGDHSLA IGTISGHARHCPDLCVWVDAHADINTPLTTSSGNLHGQPVSFLLRELQDKVPQLPGFSWIK PCISSASIVYIGLRDVPPEHFILKNYDIQYFSMRDIDRLGIQKVMERTFDLLIGKRQRPIHLSF DIDAFDPTLAPATGTPVVGGLTYREGMYIAEEIHNTGLLSALDLVEVNPQLATSEEEAKTTAN LAVDVIASSFGQTREGGHIVYDQLPTPSSPDESENQARVRI
Clonality	Polyclonal
NCBI	AAH29050.1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



ARG2 MaxPab polyclonal antibody. Western Blot analysis of ARG2 expression in COLO 320 HSR.



Immunoperoxidase of purified MaxPab antibody to ARG2 on formalin-fixed paraffin-embedded human prostate. [antibody concentration 3 ug/ml].



Western Blot analysis of ARG2 expression in transfected 293T cell line by ARG2 MaxPab polyclonal antibody. Lane 1: ARG2 transfected lysate(38.94 KDa). Lane 2: Non-transfected lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)