

Product Datasheet

Goat Anti-Rabbit IgG (H+L) (HRP) (orb229657)

Description	Goat polyclonal to Goat Anti-Rabbit IgG (H+L) (HRP)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:5000-10000
Expiration Date	12 months from date of receipt.