

Product Datasheet

CHX10 Rabbit monoclonal antibody (orb2297195)

Description	CHX10 Rabbit monoclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthesized peptide derived from human CHX10
Form/Appearance	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	50kDa
Uniprot ID	P58304
Dilution Range	WB 1:500~1:2000 FC 1:200
Expiration Date	12 months from date of receipt.