

## Product Datasheet

### FABPI Rabbit monoclonal antibody (orb2297197)

<b>Description</b>	FABPI Rabbit monoclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A synthesized peptide derived from human FABPI
<b>Form/Appearance</b>	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	15kDa
<b>Uniprot ID</b>	<a href="#">P12104</a>
<b>Dilution Range</b>	WB 1:500~1:2000
<b>Expiration Date</b>	12 months from date of receipt.