


Product Datasheet

HA-Tag antibody (orb229733)

Description	Rabbit polyclonal antibody to HA-Tag
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:3000-10000, ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of HA fusion prote...