

Product Datasheet

C4A Rabbit monoclonal antibody (orb2297449)

Description	C4A Rabbit monoclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	A synthesized peptide derived from human C4A
Form/Appearance	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	193kDa
Uniprot ID	P0COL4
Dilution Range	WB 1:500~1:2000 IHC 1:50~1:200 FC 1:50
Expiration Date	12 months from date of receipt.