



Product Datasheet

HSF1 Antibody (orb229806)


Description	Mouse monoclonal antibody to HSF1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC
Immunogen	Purified recombinant fragment of human HSF1 (AA: 256-359) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	57.3kDa
Entrez	3297
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



ELISA Black line:
Control Antigen (100
n...



Immunofluorescence
analysis of A431
cell...



Immunohistochemical
analysis of formalin...