



Product Datasheet

PER3 Antibody (orb229809)


Description	Mouse monoclonal antibody to PER3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human PER3 (AA: 723-954) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	132kDa
Entrez	8863
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



ELISA Black line: Control Antigen (100 n...



Western blot analysis of human recombina...



Western blot analysis of HEK293 using PE...