



Product Datasheet

PER3 Antibody (orb229810)


| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to PER3 |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC, ICC, IHC |
| Immunogen | Purified recombinant fragment of human PER3 (AA: 723-954) expressed in E. Coli. |
| Preservatives | Purified antibody in PBS with 0.05% sodium azide. |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | Mouse IgG2b |
| Clonality | Monoclonal |
| MW | 132kDa |
| Entrez | 8863 |
| Dilution Range | IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293 using PE...



ELISA Black line: Control Antigen (100 n...



Flow cytometric analysis of HeLa cells u...