



Product Datasheet

RAP1GAP Antibody (orb229816)


Description	Mouse monoclonal antibody to RAP1GAP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human RAP1GAP (AA: 412-660) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	73.4kDa
Entrez	5909
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



ELISA Black line: Control Antigen (100 n...



Flow cytometric analysis of A431 cells u...



Western blot analysis of human recombina...