

## Product Datasheet

# Collagen I Rabbit monoclonal antibody (orb2298372)

<b>Description</b>	Collagen I Rabbit monoclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bovine, Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	A synthesized peptide derived from human Collagen I
<b>Form/Appearance</b>	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	139kDa
<b>Uniprot ID</b>	<b>P08123, P02452</b>
<b>Dilution Range</b>	WB 1:500~1:2000 IHC 1:50~1:200
<b>Expiration Date</b>	12 months from date of receipt.