

Product Datasheet

Sulindac EP Impurity C-d3 (orb2299739)

Description	Sulindac EP Impurity C-d3 is a deuterated compound of Sulindac EP Impurity C. Sulindac EP Impurity C has a CAS number of 32004-67-4. (E/Z)-Sulindac sulfide is a potent γ -secretase modulator (GSM) that selectively reduces A β 42 production in favor of shorter A β species. (E/Z)-Sulindac sulfide can be used for researching Alzheimer's disease [1].
Storage	-20°C
Note	For research use only
Formula	C ₂₀ H ₁₄ D ₃ FO ₂ S
Purity	98.00%
MW	343.43
CAS Number	250608-66-3
SMILES	<chem>C1C=C2/C(=C\C3C=CC(SC([2H])([2H])([2H])=CC=3)/C(C)=C(CC(O)=O)C2=CC=1F</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)