

Product Datasheet

Kynurenic Acid-d5 (orb2300135)

Description	Kynurenic Acid-d5 is a deuterated compound of Kynurenic Acid. Kynurenic Acid has a CAS number of 492-27-3. Kynurenic acid, a natural metabolite of tryptophan via the kynurenine pathway, is a broad-spectrum excitatory amino acid antagonist.
Storage	-20°C
Note	For research use only
Formula	$C_{10}H_2D_5NO_3$
Purity	98.00%
MW	194.20
CAS Number	350820-13-2
SMILES	<chem>OC=1C2=C(N=C(C(O)=O)C1[2H])C(=C(C(=C2[2H])[2H])[2H])[2H]</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)