

Product Datasheet

Pirfenidone-d5 (orb2302576)

| | |
|------------------------|--|
| Description | Pirfenidone-d5 |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{12}H_{11}NO$ |
| Purity | 98.00% |
| MW | 190.25 |
| CAS Number | 1020719-62-3 |
| SMILES | <chem>O=C1N(C=C(C)C=C1)C2=C(C(=C(C(=C2[2H]))[2H]))[2H]][2H][2H]</chem> |
| Expiration Date | 12 months from date of receipt. |