

Product Datasheet

Human IL-5 Protein (orb2302587)

Description	Human IL-5 Protein, premium grade is expressed from human 293 cells (HEK293). It contains AA Ile 20 - Ser 134 (Accession # P05113-1).
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	IL-5
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 24 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 13.1 kDa. The protein migrates as 18 kDa±3 kDa under reducing (R) condition, and 36 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.
Purity	95%
MW	13.1 kDa
Uniprot ID	P05113
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)