

Product Datasheet

50bp DNA Ladder (orb2303095)

Description	The 50bp Ladder DNA Marker is composed of 9 DNA fragments with a length of 50 bp to 500 bp, dissolved in 1×Loading Buffer, and can be used for direct electrophoresis of 5-10 µl during used, which is very convenient to use. Two optimized Loading Buffer dyes are added without blocking the brightness of DNA band fluorescence.*Strips are fine-tuned, and every clip is clear enough
Note	For research use only
MW	50bp-500bp
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)