

Product Datasheet

Mouse HGF R Protein, His Tag (orb2308917)

Catalog Number	orb2308917
Category	Proteins
Description	Mouse HGF R / c-MET Protein, His Tag (orb2308917) is expressed from human 293 cells (HEK293). It contains AA Glu 25 - Ala 931 (Accession # P16056-1).
Target	HGF R
Conjugation	Unconjugated
Tag	C-10×His
Reactivity	Mouse
Form/Appearance	Powder
Buffer/Preservatives	PBS, pH 7.4
Purity	90%
UniProt ID	P16056
MW	32.1 kDa & 70.6 kDa
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 32.1 kDa (α chain) & 70.6 kDa (β chain).The protein migrates as 40-43 kDa (α subunit) and 90-95 kDa (β subunit) when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Expression System	HEK293
Endotoxins	1.0 EU per µg

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expression Region	Glu 25 - Ala 931
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)