

Product Datasheet

Human Follistatin Protein, Fc Tag (orb2321188)

Description	Human Follistatin Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gly 30 - Asn 317 (Accession # P19883-2).
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Follistatin
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 58 kDa. As a result of glycosylation, the protein migrates as 65-80 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition, and 120-130 kDa under non-reducing (NR) condition (SDS-PAGE).
Purity	90%
MW	58 kDa
Uniprot ID	P19883
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)