

Product Datasheet

RPA2 Antibody (PE/Cy5.5) (orb2334898)

| | |
|----------------------------|---|
| Description | RPA2 Antibody (PE/Cy5.5) |
| Species/Host | Rabbit |
| Conjugation | PE/Cy5.5 |
| Tested Applications | FC, WB |
| Target | RPA2 |
| Preservatives | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Recombinant |
| Uniprot ID | P15927 |
| Dilution Range | Flow-Cyt=1:50-1:100, WB=1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |