

Product Datasheet

RPA2 Antibody (PE/Cy5) (orb2334899)

Description	RPA2 Antibody (PE/Cy5)
Species/Host	Rabbit
Conjugation	PE/Cy5
Tested Applications	FC, WB
Target	RPA2
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
Uniprot ID	P15927
Dilution Range	Flow-Cyt=1:50-1:100, WB=1:500-1:2000
Expiration Date	12 months from date of receipt.