


Product Datasheet

Ferd3l antibody (orb233623)

Description	Goat polyclonal antibody to Ferd3l
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Ferd3l (aa56-68)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Frozen sections of Mouse embryo E13.5 shows staining of the floor plate of the developing spinal cord. Recommended concentration, 1-3 µg/ml.
Protein Sequence	QFDERYQEVEGDE
Clonality	Polyclonal
MW	19.5
Dilution Range	ELISA: 1:128000, IHC-P: 1-3µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Mouse emb...