



Product Datasheet

LDHB antibody (orb233628)

Description	Goat polyclonal antibody to LDHB
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	LDHB
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Reported variants represent identical proteinWB: Approx 35kDa band observed in Human, Mouse and Rat Heart lysates. Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	NARGLTSVINQKLK
Clonality	Polyclonal
MW	36.6
Entrez	3945
Dilution Range	ELISA: 1:16000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Pig Heart lysat...



Western blot analysis of Human, Mouse an...