


Product Datasheet

VRK2 antibody (orb233648)

Description	Goat polyclonal antibody to VRK2
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	VRK2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in nuclear lysates of cell line K562. Recommended concentration: 1-3 µg/ml.
Protein Sequence	PNSQKVDSQKAATK
Clonality	Polyclonal
MW	44.4; 58.1; 55.4; 45.0
Entrez	7444
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 nuclear lysates.