


Product Datasheet

FCRL1 antibody (orb233660)

Description	Goat polyclonal antibody to FCRL1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FCRL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize isoform 1WB: Approx 55+48kDa bands observed in Human Tonsil and Thymus lysates. The observed molecular weights correspond to the glycosylated forms. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	SRLRKANITDVD
Clonality	Polyclonal
MW	46.9; 40.4; 46.9
Entrez	115350
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Tonsil (L...