



Product Datasheet

TNNT2 antibody (orb233666)

Description	Goat polyclonal antibody to TNNT2
Species/Host	Goat
Reactivity	Bovine, Canine, Feline, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	cardiac troponin T (aa201-213)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize all reported isoformsWB: Approx 40kDa band observed in Human, Mouse and Pig Heart lysates
Protein Sequence	TERKSGKRQTERE
Clonality	Polyclonal
MW	35.6
Entrez	7139
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human and Mouse...



Western blot analysis of Pig Heart lysat...