

Product Datasheet

RPS6KA5 Antibody (PE/Cy5.5) (orb2337799)

Description	RPS6KA5 Antibody (PE/Cy5.5)
Species/Host	Mouse
Conjugation	PE/Cy5.5
Tested Applications	WB
Target	RPS6KA5
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1/Kappa
Clonality	Recombinant
Uniprot ID	O75582
Dilution Range	WB=1:500-1:1000
Expiration Date	12 months from date of receipt.