
Product Datasheet

PKC beta 1/PRKCB Antibody (orb234365)

Description

PKC beta 1/PRKCB Antibody


Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Western blot analysis using Anti-PKC bet...

Immunogen

E.coli-derived human PKC beta 1 recombinant protein (Position: E542-V671). Human PKC beta 1 shares 100% amino acid (aa) sequence identity with both mouse and rat PKC beta 1.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for PKC beta 1 is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

76869 MW

Uniprot ID
[P05771](#)
Expiration Date

12 months from date of receipt.