

Product Datasheet

GGT5 HC peptide (orb234491)

Description	The peptide is used to anti-GGT5 HC antibody reactivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GGT5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234860 (GGT5 HC antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P36269, Q9Z2A9, Q9QWE9
Entrez	2687, 23887
Expiration Date	12 months from date of receipt.